

ProtoCOL 3

Chromogenic ID Batch Set Up - Quick Guide

Start ProtoCOL 3 software

1

Enter



➔ Log on to ProtoCOL 3

2

Position plate

➔ Insert the plate holder, ensuring the plate holder which is recommended by the software is used

➔ Place plate onto plate holder

3

Create batch

Click



➔ Click New Batch

4

Plate configuration

Click



Defaults are Circular plate 90mm

➔ Change if necessary

5

Choose batch type

Select



➔ Select Chromogenic Media

6

Enter a Plate ID

Enter



➔ Enter a plate ID found on Measure tab. Autoincrement default is "On", to change this click on the Configuration tab

Click



Save the batch

7

Enter → Name the Batch

Click  →

Select Media Supplier and Type

8




→ Select Supplier


→ Select Agar Type

Test Measure then Measure Plate

9



Test measure to ensure user is happy with results then measure the plate



CONTACT SYNBIOSIS:

EUROPE:
 Tel: +44 (0)1223 727125
 Email: support@synbiosis.com
www.synbiosis.com

USA:
 Tel: 800 686 4451/301 662 2863
 Email: support@synbiosis.com
 P3CIDBSU1216