

ProtoCOL 3

Colour Classification Batch Set Up - Quick Guide

Start ProtoCOL 3 software

1

Enter



➔ Log on to ProtoCOL 3

2

Position plate

- ➔ Insert the plate holder, ensuring the plate holder that gives the best contrast between the colony and background is used
- ➔ Place plate onto plate holder

Capture image

3

Click



➔ Select the image tab and check the live box

Slide



➔ Adjust the exposure by increasing/ decreasing the sliding scale

Click



➔ Capture image

Create batch

4

Click



➔ Click New Batch

Plate configuration

5

Select



➔ Select Plate Configuration to set plate size, sample volume and threshold. Defaults are Circular plate 90mm, Sample volume 1ml. Change if necessary

Choose an application

6

Select



➔ Select Pour Plate

Colour classification


7

Select



➔ Select Colour Classification tab

Select



➔ Automatic Classification default is 1 colour. If a total count is required keep 1 colour selected. If differentiation by colour is required, select how many different colours are present on the plate

Classification results


8

Slide



➔ Move slider to the left to accept small particles or to right to exclude small particles

Check



➔ If any colonies are touching select the split function

Check



➔ To distinguish between two colours there may be a requirement to check Remove Fringes Around Colonies to make the colony appear as a solid colour

➔ Click Next

➔ Click Finish

Measure

9

Select



➔ Select the Measure tab

Click



➔ Test measure plate

Click



➔ Assign a plate identification template name

➔ Accept batch

Enter



➔ Name the batch

ProtoCOL 3 is now ready to make the first measurement

CONTACT SYNBIOSIS:

EUROPE:
Tel: +44 (0)1223 727125
Email: support@synbiosis.com
www.synbiosis.com

USA:
Tel: 800 686 4451/301 662 2863
Email: support@synbiosis.com
P3CCBSU1216