

Protos 3 Colour Classification - Batch set-up Quick Guide

1 Start Protos 3 software

Click



➔ Run Protos 3

2 Position plate

➔ Insert the plate holder, ensuring the plate holder which gives the best contrast between the zone and background is used

➔ Place plate onto plate holder

3 Capture image

Click



➔ Select the image tab and check the live box

Slide



➔ Adjust the exposure by increasing/ decreasing the sliding scale

Click



➔ Capture image

4 Create batch

Click



➔ Click New Batch

5 Plate configuration

Click



➔ Defaults are Circular plate 90mm
Sample volume 1mL
Change if necessary

6 Choose plate type

Select



➔ Choose a plate type


7 Colour classification

Select



➔ Select Colour Classification tab

Select



➔ Automatic Classification default is 1 colour. If a total count is required keep 1 colour selected. If differentiation by colour is required select how many different colours are present on the plate


8 Classification results

Slide




➔ Move slider to the left to accept small particles or to the right to exclude small particles

Check



➔ If any colonies are touching select the Split function

Check



➔ To distinguish between two colours there maybe a requirement to check Remove Fringes Around Colonies to make the colony appear as a solid colour

➔ Click Next

➔ Click Finish

9 Measure

Select



➔ Select the Measure tab

Click



➔ Click Test Measure Plate

Click



➔ Assign a plate identification template name

➔ Accept batch

Enter

Designing Batch:

➔ Name the batch

Protos 3 is now ready to make the first measurement

CONTACT SYNBIOSIS:
EUROPE:
Email: support@synbiosis.com
www.synbiosis.com
Protos3CCQG