

# aCOLyte

Colony Counter

## Installation on Windows Vista

### Software version 3.0.1.3d

If you are about to install this software onto a Windows Vista operating system, please take note of the following instructions.

Try to install the software as per the "standard" instructions noting that where Windows XP or Windows 2000 is referred to, Windows Vista is also included.

If you see the message "The directory name is invalid", please create a folder called Temp under Windows and try again.

If you see the error message "Windows 2000 must have at least service pack 3 loaded for correct operation, so aCOLyte cannot be installed", click OK.

A text box called the Program Compatibility Assistant appears telling you that the program might not have installed correctly.

You will now be offered the choice of:-

- 1) Re-install the software using the recommended settings, OR
- 2) This program installed correctly

Choose the first option.

Begin the installation process again by choosing to allow the software to reinstall.

Reinstall as per the instructions.

#### CONTACT SYNBIOSIS:

EUROPE  
BEACON HOUSE  
NUFFIELD ROAD  
CAMBRIDGE  
CB4 1TF

Tel: +44 (0)1223 727125  
Fax: +44 (0)1223 727101  
Email: [eurosales@synbiosis.com](mailto:eurosales@synbiosis.com)  
Email: [intlsales@synbiosis.com](mailto:intlsales@synbiosis.com)

IN USA:  
5108 PEGASUS COURT, SUITE M  
FREDERICK  
MD 21704

Tel: 800 686 4451 (toll free) / 301 662 2863  
Fax: 301 631 3977  
Email: [ussales@synbiosis.com](mailto:ussales@synbiosis.com)

[www.synbiosis.com](http://www.synbiosis.com)

