

TWO STEPS FOR FAST, ACCURATE PLATE READING



CLASSIFICATION ProtoCOL 2 Batch creation

- Differentiate between colour, size and shape
- Selection of grid templates
- Movable individual counting zones
- Separation of touching colonies
- Exclusion of unwanted items such as moulds or bubbles

MEASURE Count using ProtoCOL 2

- Automated count in seconds
- Counts recorded for the entire plate and each counting sector
- Detection of organisms as small as 43µm
- Counts and images stored automatically
- Manually add or delete colonies with an audit trail to comply with GMP/GLP
- Results can be directly transferred to a LIMS system, Excel or entered into one of ProtoCOL 2's customisable reports

 **SYNBIOISIS**
A DIVISION OF THE SYNOPTICS GROUP


Synbiosis USA Headquarters:
5108 Pegasus Court Suite M
Frederick MD 21704 USA
Tel: 800-686-4451/301-662-2863
Fax: 301-631-3977
email: ussales@synbiosis.com

Synbiosis Europe and
International Headquarters:
Beacon House Nuffield Road
Cambridge CB4 1TF UK
Tel: +44 (0)1223 727125
Fax: +44 (0)1223 727101
email: sales@synbiosis.com


Website: www.synbiosis.com

A.015.10.10 All trademarks acknowledged

OPKA - OPSONOPHAGOCYtic KILLING ASSAY




 **SYNBIOISIS**
A DIVISION OF THE SYNOPTICS GROUP



OPKA OPSONOPHAGOCYTTIC KILLING ASSAY

- The *in vitro* opsonophagocytic killing assay (OPKA) is essential for developing and improving vaccines, particularly pneumococcal vaccines
- Host protection against pneumococcal disease is mainly mediated by phagocytosis
- Opsonophagocytosis is a mechanism by which the host protects against infection, with the participation of serum opsonins (antibodies and complement)
- The presence of functional antibodies leads to an effective opsonization and recovery from infection
- OPKA makes it possible to reproducibly estimate the phagocytic titre of sera from vaccinated and unvaccinated individuals



SECTOR COUNTS



ADJUSTABLE COUNTING FRAME



REPORTING



COLONY MARKERS



ProtoCOL 2 with its easy to use software, built in MS SQL server database and touch screen interface makes the counting and comparison of OPKA plates fast, secure and accurate